


InnovaMatrix® AC

Application Guide



InnovaMatrix® AC


InnovaMatrix® AC is a significant and innovative advancement in biologic dressings for acute, traumatic, and chronic wound care. Manufactured with Convatec's proprietary TriCleanse™ Process, InnovaMatrix® AC is a new category of placental-derived medical device that preserves the inherent benefits of the placenta while adding the reliability, reproducibility, and safety of a medical device.

Indications for Use*

InnovaMatrix® AC is indicated for the management of wounds including:


* See package insert for full list of indications.

- Partial and full-thickness wounds
- Pressure ulcers
- Venous ulcers
- Diabetic ulcers
- Chronic vascular ulcers
- Tunneled/undermined wounds
- Surgical wounds (donor sites/grafts, post-Mohs surgery, post-laser surgery, podiatric, wound dehiscence)
- Trauma wounds (abrasions, lacerations, and skin tears)
- Partial-thickness second degree burns
- Draining wounds



Wound Management with InnovaMatrix® AC


To use InnovaMatrix® AC effectively, a wound assessment and routine need to be performed. The following is an overview of the process, starting with inspection of the wound.



Application Process for InnovaMatrix® AC


1. Wound Bed Preparation

- a. Before applying InnovaMatrix® AC, the wound must first be debrided of all non-viable tissue or necrotic tissue. Most importantly, ensure there is no active infection.
- b. Allow for any bleeding to stop before administering InnovaMatrix® AC.
- c. Completely and thoroughly cleanse the wound.



2. Device Size Selection

- a. Measure the wound and select the appropriate size sheet of dry InnovaMatrix® AC to cover the entirety of the wound surface.
- b. Allow for any bleeding to stop before administering InnovaMatrix® AC.




Note: Cut the sheet to cover the entire wound surface and extend slightly beyond the wound margins.

InnovaMatrix® AC	Size	Coverage
IMX-15MM-01	15mm Disc	1.77 cm ² coverage
IMX-0202-01	2 cm x 2 cm	4 cm ² coverage
IMX-0404-01	4 cm x 4 cm	16 cm ² coverage
IMX-0406-01	4 cm x 6 cm	24 cm ² coverage
IMX-0505-01	5 cm x 5 cm	25 cm ² coverage


3. Device Application

- a. Apply the InnovaMatrix® AC device that is completely dry directly on the wound to ensure ease of handling. Make sure the device is entirely in contact with the wound bed and slightly extends beyond all wound margins. If multiple sheets are needed, slightly overlap the edges of the sheets or mesh with scalpel.
- b. Securely anchor InnovaMatrix® AC with the preferred fixation method (e.g., STERI-STRIP™ tissue sealant, bolsters, dissolvable clips, sutures, staples, or other appropriate fixation method) based on the type, location, and complications of the wound.
- c. Apply sterile saline on the InnovaMatrix® AC device to rehydrate the sheet.



4. Wound Dressings

- a. To protect the device from adhering to secondary dressing, apply an appropriate non-adherent primary wound dressing over InnovaMatrix® AC.
- b. Apply appropriate secondary dressing that will manage the wound exudate and keep the InnovaMatrix® AC moist and securely in place.



5. Dressing Changes

- a. To prevent damage to the newly incorporating InnovaMatrix® AC, change the primary dressing only as necessary, typically every 7 days.
- b. Change the secondary dressing as appropriate. Take care to avoid dislodging the InnovaMatrix® AC device when secondary dressing is changed.



Management Post Initial Application

1. Wound Assessment and Wound Bed Preparation for Reapplication

- a. Change all dressings every 7 days, or as needed.


NOTE: If a gel forms on the wound surface, do not attempt to forcibly remove it. Successful absorption of InnovaMatrix® AC may form a caramel-colored or off-white gel. Do not remove this gel by debridement. This gel contains ECM, which continues to replace deficient and missing ECM in the wound.

- b. As healing occurs, sections of InnovaMatrix® AC may gradually peel around the edges. Remove any loose product as needed.
- c. Gently cleanse the wound surface with sterile saline; be sure to leave any remaining ECM gel intact.
- d. Reassess the wound and record healing progression.

2. Reapplication of InnovaMatrix® AC and Dressing Changes

- a. If the wound is free of infection and necrosis but not fully epithelialized, reapply newly prepared InnovaMatrix® AC over previously absorbed application. Reapply every 7 days or as needed.
- b. Apply new dressings as described earlier in Step 5 - Dressing Changes.

★ Please refer to
Instructions for Use (IFU)
for complete instructions.





Process Flow Chart



Material Source

Production begins with a carefully selected tissue that is sourced from a highly controlled, monitored, and exclusive facility. The site is certified to ISO 13485:2016 and ISO 9001:2015 standards and compliant with the FDA's Good Manufacturing Practices in 21 CFR 820. In addition, the raw placental tissues are compliant with ISO 22442-2 standards.



Physical Processing

The placental tissue undergoes thorough bulk washing, which removes surface contaminants, residual blood, and amniotic fluid.



TriCleanse™ Process

The placental tissue is subjected to a series of chemical baths and washes that disinfect and decellularize the tissue.



Preservation

The tissue is dried.



Final Processing

The dried tissue is processed to its final configuration.




Sterilization

The tissue is terminally sterilized to an SAL 10⁻⁶.


The InnovaMatrix[®] Equation

First placental-derived medical device
for surgical and chronic wounds

Performance



Safety



+

=

InnovaMatrix[®] Platform

- Innovative
- Reliable
- Affordable



convatec

— forever caring —

Convatec Triad Life Sciences, LLC

1770 Moriah Woods Blvd. • Suite 18 • Memphis, TN 38117 • www.triadls.com